

List of Compatible Models for Digital Color Camera DP21

We announce a list of compatible models based on the results of our compatibility testing of DP21 on various PC models.

When connecting the DP21 to the PC, confirm the IEEE1394 pin number.

6 pin connector: Need **DP21-DKT**

4 pin connector: Need **DP21-LPT**

1. Current model

1-1. Desktops

Manufacturer	Model Name	IEEE1394 Connector pin number (port number)	Necessity of external IEEE1394 board	Model name of compatible IEEE1394 board	OS	Result
DELL	OptiPlex780(MT)	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3 or BUFFALO IFC-ILP4	Windows Vista SP2	Pass
					Windows XP SP3	Pass
HP	Compaq dc7900SF/CT	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3	Windows Vista SP2	Pass
					Windows XP SP3	Pass
	Compaq 8100 Elite SF/CT	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3	Windows XP SP3	Fail
	Compaq 6000ProSF/CT	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3 or BUFFALO IFC-ILP4	Windows Vista SP2	Fail
					Windows XP SP3	Fail
Compaq 8000 Elite MT/CT	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3	Windows Vista SP2	Fail	
				Windows XP SP3	Fail	

1-2. Laptops

Manufacturer	Model Name	IEEE1394 Connector pin number (port number)	Necessity of external IEEE1394 board	Model name of compatible IEEE1394 board	OS	Result
HP	ProBook 6540b Notebook PC	4 pin (1 port)	Unnecessity	—	Windows XP SP3	Pass
	Compaq 6730b/CT	4 pin (1 port)	Unnecessity	—	Windows Vista SP1	Pass
					Windows XP SP2	Pass
	Compaq EliteBook8530w	4 pin (1 port)	Unnecessity	—	Windows Vista SP2	Fail
					Windows XP SP3	Fail
Lenovo	ThinkPad R500	4 pin (1 port)	Unnecessity	—	Windows Vista SP1	Pass
	ThinkPad R400	4 pin (1 port)	Unnecessity	—	Windows Vista SP2	Pass
DELL	Latitude E5500	4 pin (1 port)	Unnecessity	—	Windows Vista SP2	Pass
					Windows XP SP3	Pass

-DP21 can not operate correctly with some external IEEE 1394 board as below regardless of computers.

RATOC REX-PCIFW2

RATOC REX-PFW4H

- All brand names and product names are registered trademarks of their respective companies.

- The results described above do not provide any guarantees on the operation. DP21 may not operate appropriately, even if a computer has the same model name, depending on computer specifications.

2. Old model

2-1. Desktops

Manufacturer	Model Name	IEEE1394 Connector pin number (port number)	Necessity of external IEEE1394 board	Model name of compatible IEEE1394 board	OS	Result
DELL	T7400	6 pin (2 ports)	Unnecessity	—	Windows Vista SP2	Pass
					Windows XP SP2	Pass
	Optiplex 360	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3 or BUFFALO IFC-ILP4	Windows XP SP2	Pass
					Windows Vista SP2	Pass
	Vostro 220s	6 pin 2 ports(option of PC) or 3 ports	Necessity	Option of PC or I/O DATA 1394-PCI3	Windows XP SP3	Pass
	HP	Compaq dc5800SF	6 pin (3 ports)	Necessity	I/O DATA 1394-PCI3 or BUFFALO IFC-ILP4	Windows Vista SP2
Windows XP SP3						Pass

2-2. Laptops

Manufacturer	Model Name	IEEE1394 Connector pin number (port number)	Necessity of external IEEE1394 board	Model name of compatible IEEE1394 board	OS	Result
DELL	Vostro 1200	6 pin (2 ports)	Necessity	SONNET FireWire / USB ExpressCard34	Windows Vista	Pass

				(FWUSB2-E34) {1}	SP2	
					Windows XP SP2	Pass
	Vostro1710	4 pin (1 ports)	Unnecessity	—	Windows Vista SP2	Fail
					Windows XP SP3	Fail

{1} SONNET FWUSB2-E34 is the PCI express card. When using Sonnet FWUSB2-E34 express card, insert the express card before starting your computer. And FWUSB2-E34 need the AC adopter separately.

-DP21 can not operate correctly with some external IEEE 1394 board as below regardless of computers.

RATOC REX-PCIFW2

RATOC REX-PFW4H

- All brand names and product names are registered trademarks of their respective companies.

- The results described above do not provide any guarantees on the operation. DP21 may not operate appropriately, even if a computer has the same model name, depending on computer specifications