

## **ULTRASONIC TECHNOLOGY**

Rapid Cutting and Reliable Vessel Sealing from a Single Multi-Functional Instrument



## Main Features

- Wiper jaw design distributes pressure evenly along the jaw length for consistent sealing.
- Excellent grasper due to higher grasping forces at the instrument tip.
- Excellent dissector with a fine tip design and a wide opening aperture.
- Rapid tissue cutting.
- Two handle types (Inline and Pistol Grip) and four working lengths (10, 20, 35 and 45 cm) to meet a variety of surgical needs.
- Lower smoke and mist generation than leading competitors' ultrasonic device for consistently clearer surgical view.
- VAR and MAX handswitches on handpiece for surgeon control.

### Generator and accessories

**SONICBEAT**

<b>USG-400</b>	Ultrasonic Generator
<b>TC-E400</b>	Energy cart
<b>TD-SB400</b>	SONICBEAT Transducer
<b>MAJ-1869</b>	Footswitch for SONICBEAT
<b>MAJ-1871</b>	Communication Cable 0.25 m
<b>MAJ-1872*</b>	Communication Cable 10 m
<b>MAJ-1873*</b>	Adapter for UHI-2/3
<b>MAJ-1876</b>	Docking Fixture

\*Required for the automatic mist and smoke evacuation function

### SONICBEAT instruments

**SONICBEAT**

<b>SB-0545PC</b>	SONICBEAT 5 mm, 45 cm, Pistol Grip
<b>SB-0535PC</b>	SONICBEAT 5 mm, 35 cm, Pistol Grip
<b>SB-0545IC</b>	SONICBEAT 5 mm, 45 cm, Inline Grip
<b>SB-0535IC</b>	SONICBEAT 5 mm, 35 cm, Inline Grip
<b>SB-0520IC</b>	SONICBEAT 5 mm, 20 cm, Inline Grip
<b>SB-0510IC</b>	SONICBEAT 5 mm, 10 cm, Inline Grip



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.