

ULTRASONIC TECHNOLOGY

Rapid Cutting and Reliable Vessel Sealing from a Single Multi-Functional Instrument



Main Features

- Wiper jaw design distributes pressure evenly along the jaw length for consistent sealing.
- Excellent grasper due to higher grasping forces at the instrument tip.
- Excellent dissector with a fine tip design and a wide opening aperture.
- Rapid tissue cutting.
- Two handle types (Inline and Pistol Grip) and four working lengths (10, 20, 35 and 45 cm) to meet a variety of surgical needs.
- Lower smoke and mist generation than leading competitors' ultrasonic device for consistently clearer surgical view.
- VAR and MAX handswitches on handpiece for surgeon control.

Generator and accessories

SONICBEAT

USG-400	Ultrasonic Generator
TC-E400	Energy cart
TD-SB400	SONICBEAT Transducer
MAJ-1869	Footswitch for SONICBEAT
MAJ-1871	Communication Cable 0.25 m
MAJ-1872*	Communication Cable 10 m
MAJ-1873*	Adapter for UHI-2/3
MAJ-1876	Docking Fixture

*Required for the automatic mist and smoke evacuation function

SONICBEAT instruments

SONICBEAT

SB-0545PC	SONICBEAT 5 mm, 45 cm, Pistol Grip
SB-0535PC	SONICBEAT 5 mm, 35 cm, Pistol Grip
SB-0545IC	SONICBEAT 5 mm, 45 cm, Inline Grip
SB-0535IC	SONICBEAT 5 mm, 35 cm, Inline Grip
SB-0520IC	SONICBEAT 5 mm, 20 cm, Inline Grip
SB-0510IC	SONICBEAT 5 mm, 10 cm, Inline Grip



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.