

FLAT BLADED MYRINGOTOMY KNIVES AND BLADES




Angled Myringotomy Blade
70130780



Myringotomy Arrow (Juvenile) Blade
70130781




Myringotomy Lance (Juvenile) Blade
70130933




Myringotomy Spear (Juvenile) Knife
with Green Handle
70130791



Myringotomy Lance (Juvenile) Knife
with Green Handle
70130790



Tympaanoplasty Angled Blade
70130925



Tympaanoplasty Straight Blade
70130926




Stainless Steel Handle for Flat Blades
130884

ROUND SHAFT MYRINGOTOMY BLADES




Spear Micro-Point Straight Blade
70130793



Spear Micro-Point Angled Blade
70130794


ROUND BLADES WITH ADDITIONAL TAPERED SHAFTS




Tapered Series – Straight Edge Angled Up Cut
70130927




Tapered Series – Straight Edge Angled Down Cut
70130928



Tapered Series – Straight Edge Straight Blade
70130929




Tapered Series – Spear Tip Angled Blade
70130930



Tapered Series – Spear Tip Straight Blade
70130931




Stainless Steel Handle for Tapered Shaft Blades
70130932




Stainless Steel Handle for Round Shaft Blades
70130934


TRADITIONAL DISPOSABLE KNIVES WITH HANDLES



Sickle Blade Flat Handle
130715




Lancet Blade Flat Handle
130716



Sickle Blade Round Handle
130915



Lancet Blade Round Handle
130916



Arrow Blade Round Handle
70130717

OLYMPUS

For more information call 800.262.3540
For order inquiries call 800.773.4301
or visit us at www.gyrusacmi.com

ENT BUSINESS GROUP:

Gyrus ENT, L.L.C.
2925 Appling Road
Bartlett, TN 38133
TEL 1-901-373-0200
FAX 1-800-757-2942

US Headquarters:

Gyrus ACMI, Inc.
136 Turnpike Road
Southborough, MA 01772
TEL 1-508-804-2600

Sales Locations:

USA

Gyrus ACMI L.P.
6655 Wedgwood Rd, Suite 160
Maple Grove, MN 55311
TEL 1-763-416-3000

Canada

Olympus Canada, Inc.
151 Telson Road
Markham, Ontario L3R 1E7
Canada
TEL 1-800-387-0437

Latin America

Olympus Latin America
5301 Blue Lagoon Drive, Suite 290
Miami, FL 33126, United States of America
TEL 1-305-266-2352

UK

KeyMed Ltd
KeyMed House
Stock Road
Southend-on-Sea, Essex SS2 5QH
United Kingdom
TEL 44-0-1702-616333

Europe

Olympus Medical Systems Europa GmbH
Wendenstr. 14-18
20097 Hamburg
Germany
TEL 49-40-23773-0

Australia

Olympus Australia Pty Ltd.
31 Gilby Road
Mount Waverley
Victoria 3149
Australia
TEL 61-1300-132-992

Hong Kong

Olympus Hong Kong and China Ltd.
L43, Office Tower, Langham Place
8 Argyle Street, Mongkok
Kowloon, Hong Kong
TEL 85-2-2481-7812

Korea

Olympus Korea Co., Ltd.
4F, Gyeongam Bldg.
Kangnam-Gu, Seoul, 135-090
Korea
TEL 82-2-6207-3200

Singapore

Olympus Singapore PTE Ltd.
491B River Valley Road #12-01/04
Valley Point Office Tower
Singapore 248373
TEL 65-6834-0010

China

Olympus Sales & Service Co., Ltd.
R1202, NCI Tower, A12 Jianguomenwai Ave.
Chaoyang District
Beijing, 100022, China
TEL 010-65693535

OLYMPUS



MYRINGOTOMY VENT TUBE GUIDE

- FLUOROPLASTIC • PACIFIC® • SILICONE
- ULTRASIL™ • MICRON® TITANIUM
- STAINLESS STEEL • HYDROXYLAPATITE
- POLYETHYLENE

© 2009 Gyrus ENT, L.L.C. All rights reserved.

®/™ are trademarks of Gyrus Group, PLC or its subsidiaries in the U.S. and/or other countries of the world. Hydromer is a registered trademark of Hydromer Inc.

REV AA 70600977



ONLY OLYMPUS CAN OFFER GYRUS ENT'S COMPLETE LINE OF

TUBE+TONSIL™
PEDIATRIC
ENT SOLUTIONS

GYRUS ENT
Arthrocare
Medtronic

Tubes	Tonsil	
	-ectomy	-otomy
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

VENTILATION TUBE SOLUTIONS

For nearly 50 years ENT surgeons have placed their trust in our Armstrong, DuraVent™, and Modified T-Tubes, all featuring our Pacific® and UltraSil™ treatments.

TONSILLECTOMY SOLUTIONS

The TA² disposable PlasmaKnife™ tip precisely excises tonsils utilizing low RF temperature and the smallest thermal margins¹ available to help lessen pain and accelerate healing.

TONSILLOTOMY SOLUTIONS

PK diego® provides hemostasis and control using a choice of our patented PK electrosurgery powered dissection blades, or our traditional powered blades.

¹data on file

Note
Unless otherwise indicated, the left-hand side of the image indicates medial side placement.

ARMSTRONG BEVELED GROMMETS

Material	Cat. No.	ID	Description
Fluoroplastic	140242	1.14 mm	
Fluoroplastic	*145700	1.14 mm	
Fluoroplastic	**70143061	1.14 mm	Blue
Fluoroplastic	*70145061	1.14 mm	Blue
Fluoroplastic	**70145761	1.14 mm	Blue
Fluro-PacifiC	70140267	1.14 mm	
Fluro-PacifiC	**70143067	1.14 mm	
Fluro-PacifiC	*70146067	1.14 mm	
Fluro-PacifiC	70140395	1.14 mm	Arm. V Style, Blue
Fluro-PacifiC	**70140396	1.14 mm	Arm. V Style, Blue
Fluro-PacifiC	*70140397	1.14 mm	Arm. V Style, Blue
Silicone	240050	1.14 mm	w/Tab, Blue
Ultrasil	70241050	1.14 mm	w/Tab, Blue
Ultrasil	**70241350	1.14 mm	w/Tab, Blue
Ultrasil	**70245755	1.14 mm	w/Tab, Blue

Designed by B.W. Armstrong, M.D., Charlotte, NC.

ARMSTRONG BEVELED PLAIN END VENTILATION TUBES

Material	Cat. No.	ID	Length
Fluoroplastic	140244	1.14 mm	7.6 mm
Fluoroplastic	*145706	1.14 mm	7.6 mm
Fluro-PacifiC	70140345	1.14 mm	7.6 mm
Fluro-PacifiC	**70140346	1.14 mm	7.6 mm
Fluro-PacifiC	*70140347	1.14 mm	7.6 mm
Fluro-PacifiC	70140398	1.14 mm	Arm V, 4.5 mm
Fluro-PacifiC	*70140399	1.14 mm	Arm V, 4.5 mm
Fluro-PacifiC	**70140400	1.14 mm	Arm V, 4.5 mm
Fluro-PacifiC	70140401	1.14 mm	Arm V, 7.5 mm
Fluro-PacifiC	*70140402	1.14 mm	Arm V, 7.5 mm
Fluro-PacifiC	**70140403	1.14 mm	Arm V, 7.5 mm
Silicone	240052	1.14 mm	7.6 mm Blue
Ultrasil	70241052	1.14 mm	7.6 mm Blue

Designed by B.W. Armstrong, M.D., Charlotte, NC.

BAXTER BEVELED BOBBIN VENTILATION TUBE

Material	Cat. No.	ID	Description
Fluro-PacifiC	70140271	.97 mm	Blue
Fluro-PacifiC	**70140273	.97 mm	Blue
Fluro-PacifiC	*70140275	.97 mm	Blue
Fluro-PacifiC	70140270	1.27 mm	Blue
Fluro-PacifiC	**70140272	1.27 mm	Blue
Fluro-PacifiC	*70140274	1.27 mm	Blue

BERGER "V" BOBBIN VENTILATION TUBE

Material	Cat. No.	ID	Description
Titanium	145282	1.27 mm	-

Designed by Eric S. Berger, M.D., Aberdeen, WA. U.S. Patent No. 4,080,171.

BUTTERFLY VENTILATION TUBE

Material	Cat. No.	ID	Length	Width
Silicone	240073	1.27 mm	5.3 mm	7.6 mm
Silicone	**70245018	1.27 mm	5.3 mm	7.6 mm
Silicone, Ultrasil	70241073	1.27 mm	5.3 mm	7.6 mm

CANALE VENTILATION TUBE

Material	Cat. No.	ID	Description
Silicone	240080	1.14 mm	Blue

Developed in conjunction with James Canale, M.D., Blytheville, AR.

COLLAR BOBBIN® VENT TUBES

Material	Cat. No.	ID	Description
Titanium	145287	1.14 mm	-
Fluoroplastic	145288	.75 mm	-
Fluoroplastic	145289	1.14 mm	-
Titanium	145297	1.27 mm	-
Fluoroplastic	145298	1.27 mm	-

Designed for Bruce Black, M.D., Brisbane, Australia.

COLLAR BUTTON VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	145214	1.27 mm	White
Fluoroplastic	**145004	1.27 mm	White
Fluoroplastic	*145705	1.27 mm	White
Fluoroplastic	**70145062	1.27 mm	Blue
Fluoroplastic	*70145060	1.27 mm	Blue
Fluoroplastic	145217	1.27 mm	White w/Wire
Fluro-PacifiC	70145314	1.27 mm	White
Fluro-PacifiC	**70145316	1.27 mm	White
Fluro-PacifiC	*70145317	1.27 mm	White
Silicone	240014	1.27 mm	Blue
Silicone, Ultrasil	70241014	1.27 mm	Blue
Silicone, Ultrasil	**70241314	1.27 mm	Blue
Silicone, Ultrasil	*70241060	1.27 mm	Blue

DONALDSON VENTILATION TUBES

Material	Cat. No.	ID	Description
Silicone	240015	1.14 mm	Blue
Soft Silicone	240060	1.14 mm	Blue
Silicone	*245702	1.14 mm	Blue
Silicone	240016	1.14 mm	w/Tab, Blue
Silicone	**70245012	1.14 mm	w/Tab, Blue
Silicone	*245707	1.14 mm	w/Tab, Blue
Fluoroplastic	146043	1.14 mm	
Ultrasil	70241015	1.14 mm	

Designed by James A. Donaldson, M.D., Seattle, WA.

DURAVENT™ TUBES

Material	Cat. No.	ID	Description
Silicone	240075	1.27 mm	1.37 mm IF, Blue
Silicone	*245775	1.27 mm	1.37 mm IF, Blue
Ultrasil	70241075	1.27 mm	1.37 mm IF, Blue
Ultrasil	**70241375	1.27 mm	1.37 mm IF, Blue

Designed by Bruce Black, M.D., Brisbane, Australia.

DUTCHER VENTILATION TUBE

Material	Cat. No.	ID	Description
Silicone, Ultrasil	70240054	1.27 mm	4.32 mm Medial Flange Width, w/Tab, Blue

Designed for Paul Dutcher, M.D., Rochester, NY.

FEUERSTEIN SPLIT TUBE

Material	Cat. No.	ID	Length
Fluoroplastic	140248	1.02 mm	9 mm

Designed for Sidney S. Feuerstein, M.D., New York, NY.

JAHN HYDROXYLVENT® TUBES

Material	Cat. No.	ID	Description
Hydroxylapatite	140906	1.00 mm	9.5 mm Long
Hydroxylapatite	140926	1.27 mm	9.5 mm Long
Hydroxylapatite	140927	1.00 mm	9.5 mm Long w/Flp. Liner

Packaged individually.
Designed for Anthony Jahn, MD, New York, NY.

LINDEMAN-SILVERSTEIN ARROW VENTILATION TUBE

Material	Cat. No.	ID	Description
Fluoroplastic	140240	1.14 mm	-

Designed for Roger Lindeman, M.D., Seattle, WA, and Herbert Silverstein, M.D., Sarasota, FL.

MICRON™ TITANIUM BOBBIN VENTILATION TUBES

Material	Cat. No.	ID	Description
Titanium/Fluro	145220	1.00 mm	w/Flange Holes
Titanium/Fluro	145271	1.27 mm	
Titanium/Fluro	145272	1.52 mm	
Titanium	145221	0.76 mm	Infant
Titanium	145215	1.00 mm	w/Flange Holes
Titanium	145219	1.00 mm	"Unibobbin" (FH on 1 side)
Titanium	145280	1.00 mm	
Titanium	145279	1.14 mm	
Titanium	145273	1.27 mm	w/Spoon Flange
Titanium	145276	1.27 mm	+ Hydromer™ Coat
Titanium	*70145776	1.27 mm	+ Hydromer™ Coat
Titanium	*145781	1.27 mm	
Titanium	145281	1.27 mm	
Blue Titanium	70145381	1.27 mm	Blue
Blue Titanium	70145881	1.27 mm	Blue, 12/box

MINI-REVERSIBLE VENTILATION TUBE

Material	Cat. No.	ID	Description
Silicone	240042	0.89 mm	Reversible

Designed for Burton J. Cohen, M.D., Louisville, KY.

MODIFIED T-TUBES

Material	Cat. No.	ID	Length	Description
Silicone	240072	1.32 mm	4.80 mm	Blue
Silicone	**70245017	1.32 mm	4.80 mm	Blue
Silicone	*245704	1.32 mm	4.80 mm	Blue
Soft Silicone	240074	1.32 mm	4.80 mm	Blue
Ultrasil	70241072	1.32 mm	4.80 mm	Blue

MORETZ TYPE VENTILATION TUBES

Material	Cat. No.	ID	Description
Titanium	145265	1.27 mm	Angled Tab
Polyethylene	70145498	1.27 mm	Blue
Polyethylene	**70145498	1.27 mm	Blue
Polyethylene	*70145798	1.27 mm	Blue
Fluro-PacifiC	70142043	1.27 mm	

PER-LEE VENTILATION TUBE

Material	Cat. No.	ID	Description
Silicone	240026	1.57 mm	60°, 15 mm Long, Blue

Designed for John H. Per-Lee, M.D., Emory Univ., Atlanta, GA.

PAPARELLA-TYPE VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	140252	1.02 mm	w/Tab
Fluoroplastic	140254	1.27 mm	w/Tab
Silicone	240043	1.02 mm	w/oTab, Blue
Silicone	240044	1.02 mm	w/Tab, Blue
Silicone	**70245016	1.02 mm	w/Tab, Blue
Silicone	*245701	1.02 mm	w/Tab, Blue
Silicone	70140280	1.14 mm	w/Tab, Blue
Silicone	240046	1.27 mm	w/Tab, Blue
Silicone	**70243046	1.27 mm	w/Tab, Blue
Silicone	*70245746	1.27 mm	w/Tab, Blue
Silicone	240048	2.03 mm	w/Tab, Blue
Soft Silicone	240058	1.02 mm	w/Tab, Blue
Soft Sil-Ultrasil	70140281	1.14 mm	w/Tab, Blue
Soft Sil-Ultrasil	**70241380	1.14 mm	w/Tab, Blue
Soft Sil-Ultrasil	*70245780	1.14 mm	w/Tab, Blue
Silicone, Ultrasil	70241044	1.02 mm	w/Tab, Blue
Silicone, Ultrasil	**70241344	1.02 mm	w/Tab, Blue
Silicone, Ultrasil	70240280	1.14 mm	w/Tab, Blue
Silicone, Ultrasil	*70240281	1.14 mm	w/Tab, Blue
Silicone, Ultrasil	**70240282	1.14 mm	w/Tab, Blue
Silicone, Ultrasil	*145015	1.14 mm	w/Tab, Blue

Designed for Michael Paparella, M.D., Minneapolis, MN.

POPE BEVELED GROMMET VENTILATION TUBES

Material	Cat. No.	ID	Description
Polyethylene	145250	1.14 mm	Blue
Polyethylene	**70143250	1.14 mm	Blue
Polyethylene	*70145750	1.14 mm	Blue
Fluro-PacifiC	70142047	1.14 mm	Blue

Designed for T.H. Pope, Indianapolis, IN.

REUTER BOBBIN VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	145212	1.00 mm	1 mm IF, w/o Holes
Fluoroplastic	*145709	1.00 mm	1 mm IF, w/o Holes
Fluoroplastic	145213	1.14 mm	1 mm IF, w/o Holes
Fluoroplastic	**145019	1.14 mm	1 mm IF, w/o Holes
Fluoroplastic	145224	1.14 mm	1 mm IF, w/Wire
Fluoroplastic	145225	1.14 mm	1 mm IF
Fluoroplastic	*145710	1.14 mm	1 mm IF
Fluro-PacifiC	145227	1.27 mm	0.5 mm IF
Fluro-PacifiC	70145319	1.00 mm	1 mm IF w/o Holes
Fluro-PacifiC	**70145320	1.00 mm	1 mm IF w/o Holes
Fluro-PacifiC	*70145321	1.00 mm	1 mm IF
Fluro-PacifiC	70145326	1.14 mm	1 mm IF
Fluro-PacifiC	**70145327	1.14 mm	1 mm IF
Fluro-PacifiC	*70145328	1.14 mm	1 mm IF
Stainless Steel	145210	1.00 mm	1 mm IF
Stainless Steel	*145701	1.00 mm	1 mm IF
Stainless Steel	145211	1.00 mm	1 mm IF w/Wire
Stainless Steel	145216	1.00 mm	Satin, 1 mm IF w/Wire
Stainless Steel	145218	1.00 mm	Satin, 1 mm IF
Stainless Steel	145208	1.14 mm	1 mm IF
Stainless Steel	145209	1.27 mm	1 mm IF
Titanium	145215	1.00 mm	1 mm IF Micron Finish

Designed for S. Harold Reuter, M.D., Houston, TX.

RICHARDS T-TUBES

Material	Cat. No.	ID	Length	Width
Silicone	240071	1.14 mm	12 mm	9.8 mm
Silicone	**70245014	1.14 mm	12 mm	9.8 mm
Soft Silicone	240059	1.14 mm	12 mm	9.8 mm
Silicone	*245706	1.14 mm	12 mm	9.8 mm
Ultrasil	70241071	1.14 mm	12 mm	9.8 mm

ROCK PEDIATRIC VENTILATION TUBE

Material	Cat. No.	ID	Description
Fluoroplastic	140086	1.02 mm	0.89 mm IF

RUBE VENTILATION TUBE

Material	Cat. No.	ID	Length	Width	Description
Silicone	70140260	1.14 mm	4.5 mm	6.6 mm	w/Tab

SHAH VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	142230	0.76 mm	w/o Wire
Fluoroplastic	**145020	0.76 mm	w/o Wire
Fluoroplastic	146044	1.14 mm	w/o Wire
Fluoroplastic	**145001	1.14 mm	w/o Wire
Fluoroplastic	140005	1.14 mm	w/Wire
Fluoroplastic	**145015	1.14 mm	w/Wire

Designed for Mr. N. Shah, FRCS, DIO, London, UK.

SHEA PARASOL VENTILATION TUBES

Material	Cat. No.	ID	Description
Silicone	240035	1.02 mm	w/Tab, Small Flange
Silicone	240036	1.02 mm	w/Tab
Silicone	*245703	1.02 mm	w/Tab
Silicone	240038	1.52 mm	w/Tab
Ultrasil	70241036	1.02 mm	w/Tab

Designed for John J. Shea, M.D., Memphis, TN. U.S. Patent No. 3,871,380.

SHEEHY-TYPE COLLAR BUTTON

Material	Cat. No.	ID	Description
Fluoroplastic	70145971	1.27 mm	
Fluoroplastic	**70143071	1.27 mm	
Fluoroplastic	*70145771	1.27 mm	
Fluro-PacifiC	70141080	1.27 mm	White
Fluro-PacifiC	**70143075	1.27 mm	White
Fluro-PacifiC	*70145780	1.27 mm	White
Silicone	70245781	1.27 mm	Blue
Silicone	*70245782	1.27 mm	Blue

SHEPARD GROMMET VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	140007	1.02 mm	